

DingaVision Deployment Guide

Concordia Chemistry Robotics Project

Overview

This guide provides the steps to deploy and run the DingaVision web application on the Dell server. The application enables live camera input or file upload, generates a caption using the BLIP model, and allows for user-interactive interpretation of the image content. We will replace this with a more advanced model like Gemma 3 or Deepseek.

Requirements

- Windows machine (Dell preferred)
- Python 3.10+
- Access to GitHub repository

Installation Steps

1. Create Virtual Environment

```
python 3.13 -m venv .venv
```

2. Activate the Virtual Environment

- On Windows CMD:

```
.venv\Scripts\activate
```

3. Install Dependencies

```
python 3.31 -m pip install -r requirements.txt
```

4. Requirements File

The `requirements.txt` file lists all the Python packages required to run the project. To install the dependencies, run:

```
python -m pip install -r requirements.txt
```

This will install packages such as `fastapi`, `uvicorn`, `torch`, `transformers`, and `Pillow`, among others.

To Create or Update the File

If you have manually installed packages and want to generate a fresh `requirements.txt`, run:

```
pip freeze > requirements.txt
```

This captures the current environment's package list. Be cautious: this will include all packages, including transient ones, so it is recommended to clean up unused dependencies first.

Recommended Example

Below is a recommended minimal `requirements.txt` compatible with Python 3.10–3.13:

```
fastapi
uvicorn
torch
transformers
Pillow
```

Version pinning (e.g., `torch==2.1.0`) can be added for reproducibility, but it may cause incompatibility if a different Python version is used.

Run the Server

```
uvicorn main:app --host 0.0.0.0 --port 8000
```

The application will be available at: `http://localhost:8000` or via the Dell's IP on your local network (e.g., `http://192.168.50.241:8000`).

Allow Camera Permissions

- Use Google Chrome or Firefox.
- For **local connections** (`localhost`), camera access should work by default.
- For **network connections** (e.g., `192.168.x.x`), browsers may block camera access due to HTTP insecurity.
- In Firefox, set the following flags by navigating to `about:config`:
 - `media.devices.insecure.enabled = true`
- Still blocked? Use `localhost` for initial testing or configure HTTPS (not yet required for this demo).

Directory Structure

DingaVision/

static/	
index.html	% Main frontend HTML file
main.py	% FastAPI backend script
requirements.txt	% Python dependencies
.venv/	% Virtual environment

GitHub Hosting Notes

Ensure your repository includes:

- main.py
- requirements.txt
- static/index.html

Notes on Model Caching

- First-time execution will download the AI model.
- Model will be cached under `C:\Users\<username>\.cache\huggingface\hub`.
- Warning about symlinks on Windows can be safely ignored unless disk space becomes an issue.

Troubleshooting

- **Camera doesn't appear:** Make sure only one browser tab is open and that the camera isn't locked by another app.
- **Model not found or timeout:** Ensure internet is active on first launch; check for typos in model name.
- **CORS errors:** Ensure frontend is accessing the correct server port.
- **Permission errors:** Avoid running from folders with restricted access.